



ES!5; ' A

Digital Magnetic Stirrer



LABORATORY

ES-5; A

Digital Magnetic Stirrer



Cod. ES-AG3M



- Ergonomic and lightweight digital magnetic stirrer
- Adjustable speed by rotary knob
- ABS case and plate with IP42 protection class.
- Reliable device, noiseless and easy to use
- Adjustable speed up to 1500rpm
- Speed displayed on the LED digital display
- Microprocessor control technology

Technical Features:

Dimension of working plate	Ø 135 mm
Working plate material	Plastic
Motor type	DC motor
Motor rating input	5 W
Motor rating output	3 W
Power	15 W
Voltage	200–240 VAC
Frequency	50/60 Hz
Stirring positions	1
Max. stirring quantity (H ₂ O)	3 L
Max. magnetic bar (L x Ø)	50mm
Speed range	100 to 1500 rpm
Speed display	LED
Protection class according to DIN EN60529	IP42
Dimension	260x150x80 mm
Weight	1,8 Kg
Ambient temperature	5 - 40°C
Permissible relative humidity	80%

Optional Accessories:

Cod. ES-C10005	Magnetic bars, 10x6mm, box of 5 pcs
Cod. ES-C10010	Magnetic bars, 15x8mm, box of 5 pcs
Cod. ES-C10015	Magnetic bars, 30x6mm, box of 5 pcs
Cod. ES-C10020	Magnetic bars, 40x8mm, box of 5 pcs
Cod. ES-C10025	Magnetic bars, 50x8mm, box of 5 pcs
Cod. ES-C10030	Silicon protection



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

