

ESF X-RAY FILM VIEWERS

NEGATOSCOPES



X-RAY FILM VIEWERS

- Wall model. Enamelled sheet steel
- Fluorescent cold light.
- Unbreakable white monochromatic diffuser with supporting device

NEGATOSCOPE

- Modèle mural. Acier émaillé.
- Lumière froide fluorescente.
- Diffuseur blanc monochromatique incassable.

| CODE | SIZE/DIMENSIONS |
|-------------|---------------------|
| ESF2-040/00 | 38,5 x 43,5 x 12 cm |
| ESF2-080/00 | 77 x 43,5 x 12 cm |
| ESF2-120/00 | 120 x 43 x 14 cm |
| ESF2-160/00 | 160 x 43 x 14 cm |



ESF2-045/00

X-RAY FILM VIEWER

Table model
Constructed from enamelled sheet steel
Fluorescent cold light
Unbreakable white monochromatic diffuser with supporting device

NEGATOSCOPE

Modèle de table.
Acier émaillé.
Lumière froide fluorescente.
Diffuseur blanc monochromatique incassable

ESF OPTOTYPES



ILLUMINATED OPTOTYPE

Steel structure with heat toughened thermosetting powder coating
 Manufactured with acrylic screen, fuses, green light for On/Off and Shuko plug
 Viewing area 24x62 cm , external size 29x67x12 cm , weight 6 kg
 Power supply 220 V 50/60 Hz
 Same versions available on wall mounted NON ILLUMINATED chart

ECHELLE DE VISION LUMINEUX

Bâti en acier émaillé époxy
 Ecran acrylique, interrupteur on/off, fusible et fiche Shuko.
 Surface visuelle 24x62 cm , dimensions externes 29x67x12 cm , poids 6 kg
 Alimentation 220 V 50/60 Hz
 Memes versions disponibles sur mural support NON LUMINEUX

| illuminated | | Acrylic screen |
|-------------|--------------------|----------------|
| ESF2-81 | Mixed decimal 3 mt | ESF2-811 |
| ESF2-82 | Mixed decimal 5 mt | ESF2-812 |
| ESF3-83 | Monoyer 3 mt | ESF3-813 |
| ESF4-84 | Armagnac 5 mt | ESF4-814 |
| ESF5-85 | Child 3 mt | ESF5-815 |



ESSE3 srl, Via Garibaldi 30
 14022 Castelnuovo D.B. (AT)
 Tel +39 011 99 27 706
 Fax +39 011 99 27 506
 e-mail: danielec@esse3-medical.com
 web: www.esse3-medical.com

