



**DATA SHEET** 

**Model: ESB-ORTHO**

**Electric Orthopedic Traction Bed**



**Technical parameters**

External L*W*H-----	2160*1060*(500-800)mm
Mattress Platform-----	1935*900mm
Back-rest Adjustment-----	0-75°(±10°)
Knee-rest Adjustment-----	0-35°(±10°)
Reverse Trendelenburg-----	0-12°
Trendelenburg-----	0-12°

**Technical Configuration**

Linear Motor-----	4pcs
5" Double Side Secure®Castors,Central Locking-----	4pcs
ABS Side Rail-----	1set
ABS Bed Ends with Safe Lock -----	1set
IV Pole Prevision-----	4pcs
Drainage Hook-----	2pcs
Aluminium Balkan Frame-----	1set
Central Locking Pedal-----	1pc



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3.dreamgest.com](http://www.esse3.dreamgest.com)

