



ES-ELBOW Line

Elbow-operated sink



STERILIZATION



Description

Aseptic surgical scrub sink molded in one piece without any join or sealing. With no overflow hole in order to avoid bacterial development, **ES-ELBOW** has a large wall and sides protections against splashing water. Available with a full range of accessories for an ideal configuration matching your specific need. Our team advises every customer according to the country's water quality and environment.

Material: Made of reinforced molded polyester resin fiberglass. The outer coating is made of high grade iso NPG sanitary gelcoat.

Targeted services: Operating rooms, emergency boxes, intensive care unit, dental surgery offices, cosmetic and plastic surgery offices, nephrology and dialysis units, critical care, neonatal unit.

Standard equipment

- 1/2/3 Elbow control tap including mixer
- 1 Siphon of Ø40mm
- 1 Regulation valve
- HDPE flexible connecting hoses for hot and cold water
- Stainless steel fixing frame

Optional accessories

- o Prefiltering kit 1µ with pressure protection system
- o Ultraviolet kit for water disinfection - 40W
- o Medical disinfectable glass siphon
- o Terminal filter 0.2µ all germs 62 days
- o Aluminum rail for accessories + 1 clipping frame
- o Elbow command Airless Soap dispenser
- o Electronic command Airless Soap dispenser
- o Manual paper dispenser
- o Nail brush dispenser
- o Basket for accessories



ES-UNO/C

CODE

ES-UNO/C

ES-DUO/C

ES-TRIO/C

Dimensions

L : 900 mm H : 745 mm D : 570 mm Weight : 15 Kg

L : 1500 mm H : 745 mm D : 570 mm Weight : 19 Kg

L : 2000 mm H : 745 mm D : 570 mm Weight : 27 Kg



ES-DUO/C



ES-TRIO/C



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail: danielec@esse3-medical.com
web: www.esse3-medical.com

