



# ESE-i50

Full Digital B/W Ultrasound Machine



IMAGING SYSTEMS




# ESE-i50

## Full Digital B/W Ultrasound Machine

- 12.1 inch high resolution medical monitor with 90 viewing angle
- Exclusive employ double wave form imaging technology
- Outstanding portable model, compact and streamlined design
- Double activated transducer connectors, transducer auto-freeze protection
- All-round adjusting knob, efficient and fast
- One key full screen design, hidden menu item, let your diagnosis more intuitive
- Support multiple measurement software package
- Built-in lithium battery, optional trolley and suitcase
- Support multiple data transmission interface
- Auto-generation BMP format report, support USB import and export. Non use any special software to read or print on PC
- Optional laser printer can print the generated report


*Convenient workflow, excellent imaging quality*



*Trolley and suitcase*



*High density and broadband frequency transducers*



**Convex**

**Linear**

**Transvaginal**

**Micro-convex**