



ESGD-113

Electric Dermatome



LITTLE ACCESSORIES



ESGD-113

Electric Dermatome



ESGD-113

The slim line and ergonomically designed, this dermatome produces, consistent, uniform skin grafts time after time. The smooth high speed cutting action ensures an almost effortless and vibration free procedure.

Cut thicknesses are easily selected and constantly maintained using the adjustable depth gauge.

Cutting width 73 mm
Cutting thickness 0.05 - 1.00 mm
Operating speed range 2000 - 15000 rpm
Weight 0,28 kg



The dermatome head is supplied with a 73 mm blade as standard. Other blades are optional

Optional accessories



17340 - 23mm blade

17260 - 48mm blade

17270 - 73mm blade

Accurately controlled graft widths of 73 mm, 48 mm and 23 mm are readily achieved by fitting the appropriate easy-fit sole plate.



S11-103 - Sterile blades, double face (5 per box)

Stericut precision ground, double facet, ultra sharp blades provide accurate and consistent cutting performance whilst permitting very thin skin grafting.

The ESGD-113 dermatome head enables quick and safe blade replacement.



1292334 - MICROMOTOR WITH PEDAL

The dermatome head fits quickly onto the universal micro motor which is powered by the variable speed CM2 controller. The motor can be controlled either by the footswitch or directly from the controller itself. Please specify voltage and plug configuration when ordering.

Input voltage: 230/115/100 50-60Hz

Power: 60VA

Dimensions: 181 x 101 x 119 mm

Weight: 1.79/0.12 kg controller/motor



8570 - STERILISATION CONTAINER

This purpose designed compatible delivery system will ensure your instruments are maintained in good condition at all times.

The container is 100% sterilisable by steam



ESSE3 srl, Via Garibaldi 30

14022 Castelnuovo D.B. (AT)

Tel +39 011 99 27 706

Fax +39 011 99 27 506

e-mail: danielec@esse3-medical.com

web: www.esse3.dreamgest.com

