



DATA SHEET 

Model: 34086

Hypothermic blanket



Reflects and retains over 80% of radiated body heat to keep a person warm and protected even at subzero temperature. Waterproof, reusable, wind proof, lightweight. Made of polypropylene Fralom fibre coupled to a metallized polyethylene fibre. Blanket colour may change in the time.

Technical Specifications

Size: 200 x 150 cm

CE



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

