

AMBULATORY'S FURNITURE

MOD. 300445 WALL MOUNTED DRUGS CABINET

Constructed from enamelled sheet steel. Internal: one shelf Two doors with lock and key. Painting RAL 1013 ivory epoxy powder (average cover 90 micron) by industrial robot
Dimension: 75x35x45 h. 12 kg. – 0,12 cu.m.



MOD. 300450 GLASS INSTRUMENT CABINET

Constructed from sheet steel Internal: three adjustable shelves
 One glass door with lock and key. Painting RAL 1013 ivory epoxy powder (average cover 90 micron) by industrial robot Steel feet, covered with aluminium alloy tube
Dimension: 50x40x160 h. 33 kg. – 0,32 cu.m.

MOD. 300480 GLASS INSTRUMENT CABINET

Constructed from sheet steel Internal: three adjustable shelves
 Two glass doors with lock and key. Internal: one shelf
 Two doors with lock and key. Painting RAL 1013 ivory epoxy powder (average cover 90 micron) by industrial robot. Steel feet, covered with aluminium alloy tube
Dimension: 80x40x180 h. 64 kg. – 0,58 cu.m.



MOD. 300485 GLASS INSTRUMENT CABINET

Constructed from sheet steel Internal: four adjustable shelves
 Two glass doors with lock and key. Painting RAL 1013 ivory epoxy powder (average cover 90 micron) by industrial robot
 Steel feet, covered with aluminium alloy tube
Dimension: 80x40x180 h. 74 kg. – 0,58 cu.m.

MOD. 300550 WOOD MEDICINE CABINET

Realized with chipboard covered with high resistant plastic laminated material 20 mm thick.
 Drugs cabinet with door, lock and key Internal: four adjustable shelves Twenty-two shelves swivelling internal rack . Two doors with lock and key.
 Eleven shelves for each door Stainless steel socle.
Dimension: 100x60x190 h. 170 kg. – 1,14 cu.m.





Model 500690 - ADJUSTABLE HEIGHT I.V. ROD (from 150 to 220 cm.) on CASTORS
Epoxy coated and chrome - plated steel. Two hooks fitted in the holes at the four corners of the bed.
3,5kg. - 0,03cu.m.

Model 500691 - ADJUSTABLE HEIGHT I.V. ROD (from 150 to 220 cm.)
Epoxy coated and chrome - plated steel.
Two hooks fitted in the holes at the four corners of the bed.
1,5kg. - 0,01cu.m.

Model 500699 - FIXED HEIGHT I.V. ROD (150cm.)
Epoxy coated and chrome - plated steel.
Two hooks fitted in the holes at the four corners of the bed.
1,5kg. - 0,01cu.m.



MOD. 500815 WARD SCREEN

Aluminium alloy tubular frame. Two washable plastic curtains.
Dimension: 120x170 h. 5 kg. - 0,2 cu.m

MOD. 500820 WARD SCREEN

Aluminium alloy tubular frame. Three washable plastic curtains.
Dimension: 180x170 h. 9 kg. - 0,2 cu.m

MOD. 500825 WARD SCREEN

Aluminium alloy tubular frame. Four washable plastic curtains.
Dimension: 240x170 h. 12 kg. - 0,23 cu.m

MOD. 500826 CASTORS FOR WARD SCREEN

Mounted on 50 mm. castors



MOD. 500871 EXAMINATION LAMP

Aluminium alloy chrome-plated reflector adjustable in all positions from 100 to 150 cm. Cast aluminium alloy five-spoke base with five 50 mm. castors. **3 kg. - 0,02 cu.m**

MOD. 500872 HALOGEN EXAMINATION LAMP

Aluminium alloy chrome-plated reflector adjustable in all positions from 100 to 150 cm. Cast aluminium alloy five-spoke base with five 50 mm. castors.
Light intensity lux. 20 watt **4 kg. - 0,02 cu.m**



MOD. 500987 SINGLE BOWL STAND

Aluminium alloy tubular frame. Stainless steel bowl diam. 32 cm. H= 80 cm Cast aluminium five spoke-base with five 50mm. castors **8,5 kg. - 0,07 cu.m**

MOD. 500989 DOUBLE BOWL STAND

Aluminium alloy tubular frame. Stainless steel bowl diam. 32 cm. H= 80 cm Cast aluminium five spoke-base with five 50mm. castors **12 kg. - 0,10 cu.m**



MOD. 500991 WASTE BIN It. 7

Constructed from enamelled sheet steel with pedal operated lid . Removable plastic inner container.
Dimension diam. 33 cm height 38 cm
2,5 kg. - 0,03 cu.m

MOD. 500993 WASTE BIN It. 10

Constructed from stainless steel with pedal operated lid . Removable stainless steel inner container
Dimension diam. 33 cm height 38 cm
3,5 kg. - 0,04 cu.m

MOD. 500995 WALL MOUNTED X-RAY VIEWER

Epoxy powder coated steel sheet. White monochromatic diffuser.
Homogeneous fluorescent lighting. Stainless steel film holder device.
Dimension: 120x43x12 cm 14 kg. - 0,06 cu.m

MOD. 500996 WALL MOUNTED X-RAY VIEWER

Epoxy powder coated steel sheet. White monochromatic diffuser.
Homogeneous fluorescent lighting. Stainless steel film holder device.
Dimension: 43x43x12 cm 14 kg. - 0,06 cu.m

