

ESF-120051 PATIENT TROLLEY



PATIENT TROLLEY REMOVABLE TOP

Manufactured from chrome plated steel tube,
REMOVABLE TOP up-holstered and covered with washable plastic material FLAME RETARDANT
COMPLETE WITH SIDE RAILS and I.V. STAND
Mounted on 4 casters 100 mm two of which with brakes
DIMENSIONS: 180x55x80 (H) cm
Knock down construction

CHARIOT BRANCARD PLAN DE COUCHAGE AMOVIBLE

Bâti métallique, revêtement chromé,
PLAN DE COUCHAGE AMOVIBLE avec matelas fixe en tissu vinyle lavable et traité anti-feu.
LIVRÉ AVEC BARRIÈRES LATÉRALES et TIGE PORTE SÉRUM
Quatre roulettes avec \varnothing de 100 mm, dont deux avec frein.
DIMENSIONS: 180x55x80 (H) cm
Structure démontable

ESF-120050 PATIENT TROLLEY



PATIENT TROLLEY FIXED TOP

Manufactured from chrome plated steel tube
FIXED TOP upholstered and covered with washable plastic material FLAME RETARDANT.
COMPLETE WITH SIDE RAILS and I.V. STAND
Mounted on 4 casters 100 mm two of which with brakes
DIMENSIONS: 180x55x80 (H) cm Knock down construction

CHARIOT BRANCARD PLAN DE COUCHAGE FIXE

Bâti métallique, plan de couchage fixe avec matelas fixe en tissu vinyle lavable et traité anti-feu.
LIVRÉ AVEC BARRIÈRE DE PROTECTION LATÉRALES et TIGE PORTE SÉRUM
Quatre roulettes avec diam. de 100 mm, dont deux avec frein.
DIMENSIONS: 180x55x80 (H) cm
Structure démontable



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com

