



PATIENT TROLLEYS STRETCHERS





ESH-3700

Patient Stretcher Trolley

- Overall Size : 210L x 56W x 81H cms
- Frame work made of Stainless Steel
- Two section foam padded top
- Backrest adjustable
- Trolley mounted on 15cms dia Castors-two with brakes
- Removable Stretcher top



ESH-3701

Patient Stretcher Trolley

- Overall Size : 210L x 56W X 81H cms
- Frame work of Mild Steel
- Back rest Adjustable
- Collapsible safety side railings
- Cylinder Cage
- Removable stretcher top
- Pre treated and Epoxy powder coated finish
- Wire Mesh Stretcher Platform



ESH-3702

Patient Stretcher Trolley

- Overall Size : 210L x 56W x 81H cms
- Frame work of Mild Steel
- Swing away type Safety Side Railing
- Two section foam padded top
- Backrest adjustable
- Removable Stretcher Top
- Trolley mounted on 15cms dia Castors-two with brakes
- Finish : Epoxy Powder Coated



ESH-3703

Patient Stretcher Trolley

- Overall Size : 180L x 56W x 81H cms
- Frame work made of tubular M.S. Pipe
- Swing away type Safety Side Railing
- Adjustable I.V. rod, Utility Tray, Cylinder Cage
- Removable Stretcher Top
- Foam padded Top
- Trolley mount on 15 cms Dia Castors – Two with Brakes
- Epoxy Powder Coated Finish

ESH-3704

Patient Stretcher Trolley

- Overall Size : 213L x 56W x 81H cms
- Frame work of Mild Steel mounted on 15 cms castors
- Removable canvas stretcher top
- Pre treated and Epoxy powder coated finish



ESH-3705

Emergency Patient Stretcher Trolley

- Overall Size : 210L x 65W x 65-95H cms
- Safety Side Railings
- I.V. Rod
- Cylinder cage
- X-Ray translucent top with tray for holding x-ray Cassettes
- Adjustable Height
- Adjustable Height by Hydraulic Foot Operated pump
- Trendelenburg / Rev. Trendelenburg
- 15 cms dia castors
- Pre treated & Epoxy Powder Coated



ESH-3706

Emergency Patient Stretcher Trolley, Hydraulic

- Size : 1960X700X500/850mm
- Side Support : Two Plastic Moulded Track-away Safety Side Railings
- Top : Two Sections
- Framework : Mild Steel & ABS
- Finish : Beautiful Epoxy Powder Coating
- Castors : 125mm dia
- Brakes: Central Braking System
- Adjustable Functions : Height, Backrest, Trendelenburg/Reverse Trendelenburg



ESH-3707

Folding Stretcher (Canvas)

- Overall Size : 213L x 56W x 15H cms
- MildSteel Tube construction with high quality canvas material
- Can fold width wise for easy storage





ESH-3708

Folding Stretcher (Canvas)

- Overall Size : 213L x 56W x 15H cms
- Aluminium framework with high quality canvas material
- Can fold lengthwise as well as width wise

ESH-3709

Automatic Loading Stretcher for Ambulance Car



Salient Features

- This is made of high-strength aluminium alloy thick-wall pipes. It is supplied with foamed cushion and the adjustable back making the wounded feel comfortable. This machine is very convenient; one person is enough to send the patient (with this stretcher) to the ambulance car. The mechanical foldaway legs are controlled by the handles on the both ends of the stretcher. And the legs can be folded under the stretcher automatically when loading into the vehicle. Also the stretcher can be locked by the fixing devices when on the ambulance car. This type of stretcher is characterized by its light weight, durable, safe usage and easy for sterilization.
- Technical Features:
- Dimensions :
- High Position: 190x55x90cm
- Low Position: 190x60x26cm
- Max Angle of the back: 60°
- Distance between the front and rear wheel: 108cm
- Maximum Load Bearing Capacity: 159kg

ESH-3710

Aluminium Alloy Stair Stretcher

Salient Features

- This stair stretcher is made of high-strength aluminium alloy materials. It is an ideal facility for emergency staff to transfer patient to go up and down stair. Two stationary rear wheels enable the patients to be moved over most surfaces with minimal effort. The stretcher provides the patient safety during transport process because of its 3 pieces of restraint straps with quick-release buckle. Under the front of the stretcher, it has two flexible uplift rods. The stretcher can be folded and it saves the space in ambulance car. It characterized by its being light-weighted, small-sized, easy carried, user-safe and easy for cleaning.
- Technical Features
- Dimensions: 62x52x90cm
- Folding Size: 18x54x90cm
- Maximum Load Bearing Capacity: 159kg



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.esse3-medical.com

