



EST-80/DEL

*Manual and Hydraulic
Operating Table*



OPERATING THEATRE



EST-80/DEL

Manual and Hydraulic Operating Table



Available colors:



Brown



Black



Pink



Technical Features

Standard Specification

Table top size	60cm × 190cm (adds the auxiliary cart)
Elevation	70cm ~ 100cm
Angle of back rest	+ 0° ~ +50°
Trendelenburg / Anti-Trendelenburg	+ 30° ~ - 22°
Horizontal table top rotation	360°
Max. load weight	180 kg
Operation table weight	180 kg
Hand brake device	Yes

Standard Accessories

Screen frame	1 set
Shoulder support	2 pcs
Arm board	2 pcs
Foot board	2 pcs
Hand tension bar	1 pc
Hand grips	2 pcs
Auxiliary cart	1 set
Funnel	1 set
Contamination bucket	1 set



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

