



# ESF5 Line

Hand-operated vacuum devices & accessories



OPERATING THEATRE



# ESF5-09993 & ESF5-09994

## MANUAL VACUUM DEVICES



**Cod. ESF5-09993**

• **HAND-OPERATED VACUUM DEVICE**

Manual device for all suction cups characterized by a chromed steel syringe pump connected through a PVC tube to an autoclavable polycarbonate jar of 1L fitted between a metal base and a polycarbonate lid complete with vacuum gauge. Supplied in a carrying case easy to carry.



**Cod. ESF5-09994**

• **HAND-OPERATED VACUUM DEVICE**

Manual device for all suction cups characterized by a chromed steel syringe pump connected through a PVC tube to a jar of 0,5L fitted between a metal base and a lid complete with vacuum gauge. Supplied in a carrying case easy to carry.



# ACCESSORIES



Cod. **ESF5-20404** 40 mm  $\varnothing$   
Cod. **ESF5-20405** 50 mm  $\varnothing$   
Cod. **ESF5-20406** 60 mm  $\varnothing$

- Anterior suction cup, in stainless steel.



Cod. **ESF5-30700**

- Traction handle, in stainless steel.



Cod. **ESF5-10004** - 40 mm  $\varnothing$   
Cod. **ESF5-10005** - 50 mm  $\varnothing$   
Cod. **ESF5-10006** - 60 mm  $\varnothing$

- Suction cup, in stainless steel complete w/ handle.



Cod. **ESF5-6000** - 50 mm  $\varnothing$   
Cod. **ESF5-6001** - 60 mm  $\varnothing$   
Cod. **ESF5-6002** - 70 mm  $\varnothing$

- Silicone neonatal suction cup

Available in three size. in medical grade silicone rubbers. Soft suction cup that allows an easy insertion. Blue guideline designed for an easier positioning. Light weight Polysulfone handle. Autoclavable at 134°.



# ACCESSORIES



- Cod. **ESF5-25006** 6 mm  $\varnothing$
- Cod. **ESF5-25008** 8 mm  $\varnothing$
- Cod. **ESF5-25010** 10 mm  $\varnothing$
- Cod. **ESF5-25012** 12 mm  $\varnothing$
- Cod. **ESF5-25014** 14 mm  $\varnothing$

• Suction cannula, in stainless steel.

Cod. **ESF5-25015**

• Connector

Cod. **ESF5-2500**

• SET FOR UTERINE ASPIRATION CANNULA

Made of stainless steel with suction regulation handle.  
Cannule  $\varnothing$  mm. 6/8/10/12/14



- Cod. **ESF5-7030** 6 mm  $\varnothing$ , 27 cm length
- Cod. **ESF5-7031** 8 mm  $\varnothing$ , 27 cm length

• Yankauer suction cannula with handle,  
in stainless steel.



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

