



ES-HAR900

Portable auto refractometer



OPHTHALMOLOGY / ENT



ES-HAR900

Portable auto refractometer

Features

ES-HAR900 is a portable vision screening instrument used in all age groups. The working principle of the ray through the eye refraction refractive system focusing on the sensor, the signal processing through the eyes of the spherical, cylindrical diopter and axial refractive data.

1. *One key operation for test result*
2. *One measurement detect multiple eye problems, including myopia, hyperopia, astigmatism and anisometropia.*
3. *Accurate results in less than 5 seconds.*
4. *Intelligent flashing light and sound guidance requires minimum patients' cooperation. It is especially suitable for young children.*
5. *It can be widely used in vision screening for children, adults, and seniors, even patients wearing glasses or contact lenses.*

The advanced design concept combines ergonomics and aesthetic. Easy and comfortable operation is achieved with a retractable hand strap

Specification

Sphere range: Myopia Mode: -9.00D ~ +2.50D (Interval: 0.25D)

Normal Mode: -6.00D ~ +8.00D (Interval: 0.25D)

Sphere tolerance: ±0.50D or ±20% (whichever is greater)

Cylinder range: -3.00D ~ 0.00D (Interval: 0.25D)

Cylinder tolerance: ±0.50D or ±20% (whichever is greater)

Axis: 1° ~ 180° (Interval: 1°)

Axis tolerance: ±5°

Working distance: Myopia Mode: 28cm

Normal Mode / Infant mode: 35cm

Dimensions: 200mm (length) × 80mm (width) × 270mm (height)

Weight: 0.9kg

Power: 100-220V/50-60HZ/1.5A

Charge time: ≤ 2 hours

Battery life: ≤ 4 hours

Printer: Portable thermal printer



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

