



ES-OPH-360

OPHTHALMIC UNIT

UNITE' D'OPHTALMOLOGIE



2-3 instruments

Rotating table

Rotating chair



THE CHAIR IS
SELECTABLE



ES-OPH-360

OPHTALMIC UNIT

Adopt modern Italian style in appearance designing

- Human-oriented design, complete functions makes operation be more conveniently
- The chair could up and down, rotate left and right, and with tiltable backrest

Specification:

- Unit size: 1200*1055*1850mm
- Table Size: 800*490 mm
- Table rotation / movable ranges: 90°/ 335 mm
- Swing arm oscillation angle: $\pm 35^\circ$
- Chair up and down range: 540-690mm
- Chair loading ability: 250kg ·
- Adjustable backrest range: 160°
- Lamp: LED (24V, 3W)
- Electric Voltage: AC220V/110V 50/60HZ
- Weight: 200kg

UNITE' D'OPHTALMOLOGIE

Adoptant un design italien modern et fonctionnel

- Design « human-oriented » pour une facilité d'emploi opérationnelle
- La chaise peut monter et descendre, tourner à gauche et droite, avec appui bras relevables

Spécification

- Dimensions unité 1200*1055*1850 mm
- Dimensions de la table 800*490 mm
- Rotation / Mouvements de la table 360° / 335 mm
- Oscillation des bras : $\pm 35^\circ$
- Range élévation de la chaise 540-690 mm
- Charge max admissible : 250 kg
- Réglage inclinaison dossier : 160°
- Lampe visite intégrée : LED (24V, 3W)
- Alimentation 220/110V 450/60 Hz
- Poids lourd 200 kg



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.esse3-medical.com

