



ES-LS5

Digital Slit lamp



OPHTHALMOLOGY / ENT



Main features

High Definition Integrated System. Professional Software. Joystick Handle.

*This Product Must Be Connected with Computer.

Optical System: ESLM-2 Based

Digital System:

3.0m Max. Pixels Resolution, Image Capture Card, Professional Image Software

It based on our standard slit lamp and adopts CANON 1300D digital camera with which allows to take very detailed part the patient eye.

-Its professional software system includes high definition video function (MOV format), doctors can record, play, capture image from the video via software directly any time you want during playing.

-The SLR digital Camera adopts AC power to supply the power. It does not need battery.

This function is very convenient for doctors for using.

LED LAMP Optional

Specification

Type	Parallel Galilean
Eyepiece	12.5X
Total Magnification	6X, 10X, 16X, 25X, 40X
Diameter of Visual Field	Φ 37, Φ 23, Φ 14, Φ 8.7, Φ 5.7mm
Diopter Adjustment	+5D--5D
Slit Width	0mm-14mmcontinuously adjustable
Slit Height	1mm-14mmcontinuously adjustable
Aperture Diameter	Φ 14, Φ 10, Φ 5, Φ 3, Φ 2, Φ 1, Φ 0.2mm
Slit Angle	0° -180°
Illumination Tilting	0° , 5° , 10° , 15° , 20°
Filter	Heat absorption, grey, redfree(green), cobalt blue
Input Voltage	110V/220V 60/50Hz
Input Power	60VA
Illumination Bulb	,led lamp
Fixation	Red LED



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

